




S P O R T S M E D I C I N E

Force Fiber®

High Strength Suture





Tornier Force Fiber® Suture

Composed of a unique, ultra-high-molecular-weight polyethylene fiber, Tornier Force Fiber® Suture offers the strength and flexibility needed for everything from minor cuts to major surgery, including approximation and ligation of soft tissues – especially the use of allograft tissue for orthopaedic procedures. Available in size 2 and 3-4, blue/white or black/white cobraid, these high strength sutures are available with or without attached needles, come preloaded in all Insite™ suture implants, or packaged as a suture product ready for surgery.

- Ultra-high-molecular-weight polyethylene fiber provides flexibility and unmatched strength.
- Unique braid construction has 47% better knot break strength than Fiberwire™.
- Lubricious design allows each knot to slide for smoother, more efficient performance.
- A flatter surface allows each knot to lock against the previous, providing superior knot security.

Ordering Information

CAT#	Material	Size	Length	Color	Tipped	Needle	Dist. Config.
SMS100101	UHMWPE	Size 2	36 in.	Blue/White Co-braid	NO	KHCS (K-Point)	Free Pack
SMS100201	UHMWPE	Size 2	36 in.	Black/White Co-braid	NO	KHCS (K-Point)	Free Pack
SMS100301	UHMWPE	Size 3-4	36 in.	Blue/White Co-braid	NO	Diamond Point	Free Pack
SMS100401	UHMWPE	Size 3-4	36 in.	Black/White Co-braid	NO	Diamond Point	Free Pack
SMS100501	UHMWPE	Size 2	36 in.	Blue/White Co-braid	NO	HCS (Tapered)	Free Pack
SMS100601	UHMWPE	Size 2	36 in.	Black/White Co-braid	NO	HCS (Tapered)	Free Pack
SMS100701	UHMWPE	Size 3-4	36 in.	Blue/White Co-braid	NO	HCS (Tapered)	Free Pack
SMS100801	UHMWPE	Size 3-4	36 in.	Black/White Co-braid	NO	HCS (Tapered)	Free Pack
SMS100901	UHMWPE	Size 2	36 in.	Blue/White Co-braid	NO	N/A	Free Pack
SMS101001	UHMWPE	Size 2	36 in.	Black/White Co-braid	NO	N/A	Free Pack
SMS101101	UHMWPE	Size 3-4	36 in.	Blue/White Co-braid	NO	N/A	Free Pack
SMS101201	UHMWPE	Size 3-4	36 in.	Black/White Co-braid	NO	N/A	Free Pack
SMK100103	UHMWPE	Size 2	36 in.	Blue/White Co-braid	NO	N/A	Box of Twelve
SMK100203	UHMWPE	Size 2	36 in.	Black/White Co-braid	NO	N/A	Box of Twelve
SMT100103	UHMWPE	Size 2	36 in.	Blue/White Co-braid	YES	N/A	Box of Twelve
SMT100203	UHMWPE	Size 2	36 in.	Black/White Co-braid	YES	N/A	Box of Twelve
SMS100103	UHMWPE	Size 2	36 in.	Blue/White Co-braid	NO	KHCS (K-Point)	Box of Twelve
SMS100203	UHMWPE	Size 2	36 in.	Black/White Co-braid	NO	KHCS (K-Point)	Box of Twelve
SMS100303	UHMWPE	Size 3-4	36 in.	Blue/White Co-braid	NO	Diamond Point	Box of Twelve
SMS100403	UHMWPE	Size 3-4	36 in.	Black/White Co-braid	NO	Diamond Point	Box of Twelve
SMS100503	UHMWPE	Size 2	36 in.	Blue/White Co-braid	NO	HCS (Tapered)	Box of Twelve
SMS100603	UHMWPE	Size 2	36 in.	Black/White Co-braid	NO	HCS (Tapered)	Box of Twelve
SMS100703	UHMWPE	Size 3-4	36 in.	Blue/White Co-braid	NO	HCS (Tapered)	Box of Twelve
SMS100803	UHMWPE	Size 3-4	36 in.	Black/White Co-braid	NO	HCS (Tapered)	Box of Twelve
SMS100903	UHMWPE	Size 2	36 in.	Blue/White Co-braid	NO	N/A	Box of Twelve
SMS101003	UHMWPE	Size 2	36 in.	Black/White Co-braid	NO	N/A	Box of Twelve
SMS101103	UHMWPE	Size 3-4	36 in.	Blue/White Co-braid	NO	N/A	Box of Twelve
SMS101203	UHMWPE	Size 3-4	36 in.	Black/White Co-braid	NO	N/A	Box of Twelve

Distributed by: Tornier, Inc.
Edina, MN 55435 +1 888 867 6437 +1 281 494 7900
www.tornier-us.com

Tornier is a registered trademark of Tornier, SA.

This product is covered by US patent: Force Fiber Patent Pending. Additional patent(s) pending.

P/N: BAW-0730 REV:A

TORNIER 