

piton (pronounced pee'-ton)

a metal spike that is driven into a crack or seam in rock with a hammer, which acts as an anchor to protect the climber against a fall or to aid in climbing.

The Piton Knotless Fixation Implant:

- Completely eliminates knot tying
- Direct insertion; requires no tapping, punching, or pre-drilling
- Sutures can be independently adjusted and tensioned, even after instrument removal
- Excellent for both single row or double row repairs



One complete system for multiple applications

Allowing the physician to address multiple surgical procedures, the Piton™ can be used in hard and soft bone. Single and double row repairs can now be addressed entirely without cumbersome knot tying using simple or mattress suturing technique. Consider the Piton knotless fixation implant for:

SHOULDER

Rotator Cuff Repair, Bankart Repair, SLAP Lesion Repair, Biceps Tenodesis, Acromio-Clavicular Separation Repair, Deltoid Repair, Capsular Shift or Capsulabral Reconstruction

FOOT/ANKLE

Lateral Stabilization, Medial Stabilization, Achilles Tendon Repair, Hallux Valgus Reconstruction, Mid-foot Reconstruction, Metatarsal Ligament Repair

KNEE

Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Patellar Tendon Repair, Posterior Oblique Ligament Repair, Iliotibial Band Tenodesis

HAND/WRIST

Scapholunate Ligament Reconstruction, Ulnar Collateral Ligament Reconstruction, Radial Collateral Ligament Reconstruction



S P O R T S M E D I C I N E

piton™

Knotless Fixation System



Ordering Information

PITON KNOTLESS ANCHORS (SINGLE)

CAT#	Size	Material
SMK000101	3.5 mm	NiTi
SMK000201	3.5 mm	NiTi

PITON KNOTLESS ANCHORS (BOX OF 6)

SMK000503	3.5 mm	NiTi
-----------	--------	------

SUTURE INFORMATION

Size	Type	Material	Color
#2	Force Fiber	UHMWPE	Blue/White Co-braid
#2	Force Fiber	UHMWPE	Black/White Co-braid
#2	Force Fiber	UHMWPE	Blue/White Co-Braid (3), Black/White Co-braid (3)

Distributed by: Tornier, Inc.
Edina, MN 55435 +1 888 867 6437 +1 281 494 7900
www.tornier-us.com
Tornier is a registered trademark of Tornier, SA.
This product is covered by US patent: Piton Patent Pending. Additional patent(s) pending.
P/N: BAW-0732 REVA

TORNIER



